



### MODIFIED HCG TRIGGER DOSES

1000 iu HCG - draw up and mix with 3mL of diluent - administer 0.3mL of HCG

1500 iu HCG - draw up and mix with 2mL of diluent - administer 0.3mL of HCG

2000 iu HCG - draw up and mix with 2.5mL of diluent - administer 0.5mL of HCG

2500 iu HCG - draw up and mix with 2.0mL of diluent - administer 0.5mL of HCG

3300 iu HCG - draw up and mix with 3.0mL of diluent - administer 1.0mL of HCG

5000 iu HCG - draw up and mix with 2.0mL of diluent - administer 1.0mL of HCG

6700 iu HCG - draw up and mix with 1.5mL of diluent - administer 1.0mL of HCG

7500 iu HCG - draw up and mix with 2.0mL of diluent - administer 1.5mL of HCG

10,000 iu HCG - draw up and mix with 1.5mL of diluent - administer 1.5mL of HCG